

Señores:
PATRIMONIO AUTÓNOMO CENIT CONVENIOS
Ciudad

Ref.: TEK-C- 134 -23
Fecha: 28 de agosto de 2023

Referencia: MECANISMO DE ELECCION COMPETITIVO CERRADO No. 002 -2023

De acuerdo con el informe de evaluación nos permitimos enviar los siguientes documentos:

1. Solicitud de la entidad

Anexar documento de casa fabricante, en donde se especifique, que la licencia del software debe ser de última versión y venir a nombre de la Policía Nacional, además de ser a perpetuidad.

Respuesta Innovatek

Adjunto encontrarán la confirmación del fabricante sobre el cumplimiento de los requisitos establecidos para el software. Es de anotar que en los pliegos no era específico la comunicación del fabricante.

Cómo proveedores del equipo se asegura el cumplimiento del requisito de entregar el software con última versión a nombre de policía nacional

2. Solicitud de la entidad

No se evidencia que el detector del FTIR con sulfato de trigliceno, se encuentre dopado con L-alanina deuterada DLATGS.

Respuesta Innovatek

Confirmamos que la tecnología del detector es dLaTGS, cómo se evidencia en el documento adjunto especificando que los detectores de Thermo Scientific cumplen con esta tecnología desde 2007.

3. Solicitud de la entidad

No se evidencia carta del fabricante, en donde se relacione aprobaciones regulatorias de la Comisión Europea (CE) o intertek listed (ETL) o TÜV Rheinland.



Respuesta Innovatek

Se anexa certificado de fábrica, declaración de conformidad, en el cual se puede verificar el cumplimiento de las aprobaciones solicitadas

A handwritten signature in black ink, appearing to read "Hislen Hernández Rodríguez", written over a horizontal line.

Hislen Hernández Rodríguez
Representante legal
Innovatek SAS



THERMO FISHER SCIENTIFIC END-USER SOFTWARE LICENSE AGREEMENT

Please read this document carefully. This End-User Software License Agreement (the "Agreement") is a legal contract between you (hereafter in this Agreement "you" or "Licensee") and Thermo Electron Scientific Instruments LLC ("Thermo Fisher Scientific"). By installing, breaking any seals of, and/or using the Thermo Fisher Scientific molecular spectroscopy and microanalysis instrumentation software and the published user documentation accompanying the software (the "Software"), you agree to be bound by the terms and conditions of this Agreement, including without limitation the license restrictions, warranty disclaimers and limitations of liability. Subject to payment of applicable license fees, Licensee shall own the media on which Thermo Fisher Scientific originally provides the Software to Licensee hereunder (if media is provided), or on which Licensee subsequently records Software, but Thermo Fisher Scientific shall retain ownership of all copies of Software or portions thereof embodied in such media.

1. TERMS AND CONDITIONS OF SOFTWARE LICENSE.

1A. **GRANT OF LICENSE; USE RESTRICTIONS.** In consideration for the payment by Licensee of a license fee (which in some circumstances may be a part of the price Licensee paid for the Thermo Fisher Scientific instrument purchased by Licensee) or a services fee as specified in Thermo Fisher Scientific's quotation ("Quotation"), Thermo Fisher Scientific grants to Licensee a revocable, nontransferable (except as provided below) and nonexclusive right to use the Software for the number of users and in the manner permitted by the type of license for which the corresponding fee has been paid by Licensee as described in Section 1B. below and as specified in the Quotation, solely for Licensee's internal business purposes (the "License"). The term of the License is specified in the Quotation. Under no circumstances may the Software be used in any other manner or for any other purpose. Licensee agrees that it shall not sell, lease, sublicense, distribute, modify, reverse compile, decrypt, extract, disassemble, or otherwise reverse engineer any portion of the Software. Prior to disposing of any media or instrument containing Software, Licensee will ensure that any Software contained by such media or stored in such instrument has been completely erased or otherwise destroyed.

1B. TYPES OF LICENSES AND APPLICABLE RESTRICTIONS.

Licensee shall comply with the applicable License terms below based on the type of License purchased by Licensee.

- i) Single User Online License – Single user license for two computers in connection with use of the applicable instrument. The Software may reside only on two specific workstation computers ("Designated Machines") residing at the instrument location, one connected to the instrument and one on a desktop. For this type of License, no other copies are allowed for additional workstations owned or operated by the single user or Licensee. One copy may be retained for archival purposes.
- ii) Single User Offline License – Single user license for "offline" data processing (i.e., the computer is not connected to an instrument). The Software may reside on a single Designated Machine only being used for offline data processing. Under this License no other copies are allowed for additional workstations owned or operated by the single user or Licensee. One license, one computer rules apply. One copy may be retained for archival purposes. If a Single User Online License has been obtained, an additional Single User Offline License must be obtained for offline use on a designated separate computer.
- iii) Multi User License - Multi User Software. The Multi User License permits the use of one copy of the Software program or data files on a single network server ("Designated Primary Machine") and a pre-specified limited number of program files on individual client workstations, including laptop computers ("Designated Secondary Machine(s)"), in accordance with the specific licensing terms of this Agreement. These Designated Secondary Machines are

connected to the Designated Primary Machine by cable, fiber optic connection, wire, wireless, or other means and are installed with the Software. A network does not consist of computers connected to the Internet or World Wide Web but a single area network entity within a site or physical location. The Software is "in use" on a computer when it is loaded into the temporary memory (i.e. RAM) or installed into the permanent memory (e.g., DVD, hard disk, or other data storage device) of that computer. Both the program and data files must reside on the physical Local Area Network (LAN).

- iv) Cloud License. The Software is hosted by Thermo Fisher Scientific and Licensee has access to the Software via an internet connection. Licensee's use of the Software under this License is subject to the restrictions contained on the Quotation.

2. COPY RESTRICTIONS. Licensee may make one copy of the Software for backup or archival purposes to support the right to use the License; provided, however, that in no event may any copies of Software be stored on more than one workstation, computer or server at the same time. Licensee agrees that all copies of the Software shall be owned by Thermo Fisher Scientific or its suppliers and licensors, shall be considered Software subject to the terms of this Agreement and shall include the same proprietary and copyright notices and legends as supplied by Thermo Fisher Scientific and its suppliers and licensors.

3. OWNERSHIP OF SOFTWARE; CONFIDENTIALITY. Licensee agrees that no title to the Software, or the intellectual property in any of the Software, or in any Software copy, is transferred to Licensee, and that all rights in the Software and enhancements, derivatives, modifications, and improvements thereto, are owned and reserved by Thermo Fisher Scientific and its suppliers. All rights not expressly granted to Licensee under this Agreement are reserved by Thermo Fisher Scientific. This License is not a sale of the original Software or any copy thereof. Licensee shall not alter or remove any copyright, trademark or other protective or proprietary notices from the Software or media in which it is contained.

The Software, in object and source code form, the user documentation, and any related technology, idea, algorithm, intellectual property or information contained therein, including without limitation any trade secrets related to any of the foregoing, as well as any copies and derivative works, are the confidential and proprietary information ("Confidential Information") of Thermo Fisher Scientific and its suppliers. Licensee shall protect the Confidential Information from unauthorized dissemination, using the same degree of care that the Licensee ordinarily uses with respect to its own proprietary information, but in no event with less than reasonable care. Licensee shall maintain the confidentiality of the Software indefinitely. Licensee shall not use the Confidential Information for any purpose not expressly permitted by this Agreement, and shall limit the disclosure of and access to the Confidential Information to the employees or agents of Licensee who have a need to know for purposes of this Agreement, and who are, with respect to the Confidential Information, bound in writing by confidentiality terms no less restrictive than those contained herein.

4. TRANSFER RESTRICTIONS. If Licensee transfers ownership, or otherwise disposes, of a Designated Machine, provided that the transferee agrees to accept and comply with the terms and conditions of this Agreement, Licensee may transfer the Software and all licenses and rights in Software granted under this Agreement to such transferee (provided that all Software copies are also transferred). Except as provided above, or with Thermo Fisher Scientific's prior written consent, neither this Agreement nor any rights or obligations under this Agreement, in whole or in part, shall be sublicensed, assigned or otherwise transferred and any attempt to sublicense, assign or transfer this Agreement or any rights or obligation under this Agreement shall be null and void (except as aforesaid).

5. EXPORT RESTRICTIONS. Licensee acknowledges that the Software may be subject to export controls of the U.S. government. The export controls may include, but are not limited to, those of the Export Administration Regulations of the U.S. Department of Commerce (the "EAR"), which may restrict or require licenses for the export of Software from the United States and their re-export from other countries. Licensee shall comply with the EAR and all other applicable laws, regulations, laws,

treaties, and agreements relating to the export, re-export, and import of the Software. Licensee shall not, without first obtaining the required license to do so from the appropriate U.S. government agency: (i) export or re-export any Software, or (ii) export, re-export, distribute or supply any Software to any restricted or embargoed country, person or entity. Licensee shall cooperate fully with Thermo Fisher Scientific in any official or unofficial audit or inspection related to applicable export or import control laws or regulations, and shall indemnify and hold Thermo Fisher Scientific harmless from, or in connection with, any violation of this Section by Licensee or its employees, consultants, agents, or customers.

6. **TERMINATION.** If Licensee breaches this Agreement, Thermo Fisher Scientific may, in addition to any other available remedies, terminate this Agreement and all licenses and rights granted to Licensee under this Agreement. Upon the termination of this Agreement, Licensee shall, within thirty (30) days after termination, delete or destroy the Software, all removable media and documentation containing the Software.

7. **GOVERNING LAW; DISPUTES.** This Agreement and any dispute arising hereunder shall be governed by and interpreted in accordance with the laws of the State of Wisconsin, excluding its choice of laws rules. The parties hereby agree that any dispute regarding the interpretation or validity of, or otherwise arising out of, this Agreement, or relating to the products licensed hereunder shall be subject to the exclusive jurisdiction of the Wisconsin state courts of Dane County, Wisconsin (or, if there is federal jurisdiction, the United States District Court for the State of Wisconsin), and the parties agree to submit to the personal and exclusive jurisdiction and venue of these courts. If Licensee is a state or state instrumentality, Licensee hereby expressly waives the benefits of the Eleventh Amendment to the Constitution of the United States, including without limitation its application to any claim or counterclaim brought under United States Copyright or Patent Laws.

8. **INTEGRATION.** This Agreement is the entire agreement between Licensee and Thermo Fisher Scientific relating to Software and: (i) supersedes all prior or contemporaneous oral or written communications, proposals and representations with respect to its subject matter; and (ii) prevails over any conflicting or additional terms of any quote, order, acknowledgment, or similar communication between parties during the term of this Agreement. No modification to the Agreement will be binding, unless in writing and signed by a duly authorized representative of each party.

9. **FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES.** You acknowledge that the Software has not been cleared, approved, registered or otherwise qualified (collectively, "Approval") by Thermo Fisher Scientific with any regulatory agency for use in diagnostic or therapeutic procedures, or for any other use requiring compliance with any federal or state law regulating diagnostic or therapeutic products, blood products, medical devices or any similar product (hereafter collectively referred to as "federal or state drug laws"). The Software may not be used for any purpose that would require any such Approval unless proper Approval is obtained. You agree that if you elect to use the Software for a purpose that would subject you or the Software to the jurisdiction of any federal or state drug laws, you will be solely responsible for obtaining any required Approvals and otherwise ensuring that your use of the Software complies with such laws.

10. **LIMITED WARRANTY AND DISCLAIMER OF WARRANTY.**

Thermo Fisher Scientific warrants that the Software will substantially conform to the specifications which are contained in the user documentation that accompanies the Software. Thermo Fisher Scientific's warranty and obligation shall extend for a period of ninety (90) days from the date of the original delivery of Software to Licensee by Thermo Fisher Scientific ("Warranty Period"), and is solely for the benefit of Licensee, who has no authority to assign or pass through this warranty to any other person or entity. In the event Licensee notifies Thermo Fisher Scientific in writing of non-conforming Software during the Warranty Period, Thermo Fisher Scientific's entire liability and Licensee's exclusive remedy shall be for Thermo Fisher Scientific to replace the non-conforming Software pursuant to Thermo Fisher Scientific's then current terms for return and replacement of non-conforming Software.

Except as provided in this Section, the Software is provided to Licensee on an "AS IS" basis, and Thermo Fisher Scientific makes no other warranty of any kind, express or implied, with regard to the

Software. Thermo Fisher Scientific does not warrant or represent that the Software will operate uninterrupted or error free or that all defects in the Software are correctable or will be corrected. This warranty shall not apply if Software is used other than in accordance with Thermo Fisher Scientific's written instructions as contained in the user documentation, or if any of Licensee's hardware equipment associated with the use of Software malfunctions, or if the Software non-conformity is caused by an act or omission of Licensee of a third party or an event outside of Thermo Fisher Scientific's control. Software is not designed or licensed for use in on-line control equipment in hazardous environments such as operation of nuclear facilities, aircraft navigation or control, or direct life support machines.

THE OBLIGATIONS CREATED BY THIS SECTION TO REPLACE NON-CONFORMING SOFTWARE SHALL BE THE SOLE REMEDY OF LICENSEE IN THE EVENT OF NON- CONFORMING SOFTWARE. THE FOREGOING WARRANTIES ARE IN LIEU OF, AND THERMO FISHER SCIENTIFIC DISCLAIMS, ALL OTHER WARRANTIES, EXPRESS OR IMPLIED, INCLUDING, WITHOUT LIMITATION, ANY WARRANTY OF NONINFRINGEMENT, MERCHANTABILITY OR FITNESS FOR A PARTICULAR PURPOSE. THERMO FISHER SCIENTIFIC EXPRESSLY DISCLAIMS ANY WARRANTY WHATSOEVER WITH RESPECT TO ANY THIRD PARTY PRODUCTS WHETHER PROVIDED TO LICENSEE DIRECTLY BY A THIRD PARTY VENDOR OR SUBLICENSSED OR RESOLD TO LICENSEE BY THERMO FISHER SCIENTIFIC; PROVIDED, HOWEVER, THAT THERMO FISHER SCIENTIFIC WILL ASSIGN SUCH WARRANTIES TO LICENSEE WHERE POSSIBLE.

11. LIMITATION OF LIABILITY; INDEMNIFICATION. IN NO EVENT SHALL THERMO FISHER SCIENTIFIC OR ITS LICENSORS BE LIABLE TO LICENSEE OR ANY OTHER PARTY FOR ANY SPECIAL, DIRECT, CONSEQUENTIAL, INCIDENTAL, EXEMPLARY OR INDIRECT DAMAGES OF ANY KIND (INCLUDING WITHOUT LIMITATION LOSS OF PROFITS, INFORMATION OR DATA AND PERSONAL INJURY) ARISING OUT OF OR RELATED TO THIS AGREEMENT, WHETHER OR NOT THERMO FISHER SCIENTIFIC HAS BEEN ADVISED OF THE POSSIBILITY OF SUCH LOSS, HOWEVER CAUSED AND ON ANY THEORY OF LIABILITY. THESE LIMITATIONS SHALL APPLY NOTWITHSTANDING THE FAILURE OF THE ESSENTIAL PURPOSE OF ANY LIMITED REMEDY.

In no event will Thermo Fisher Scientific be liable for any claim against Licensee by a third party for any reason, and Licensee hereby agrees to defend, indemnify and hold Thermo Fisher Scientific harmless for any claims, cost, damage, expense or liability arising out of or in connection with the installation, use and performance of the Software, whether alone or in combination with any other product or service.

12. U.S. GOVERNMENT RESTRICTED RIGHTS LEGEND. The Software provided under this Agreement are commercial computer software programs developed exclusively at private expense. Use, duplication or disclosure of the Software by the Government is subject to restrictions as set forth in FAR 52.227-19(c) or subparagraph (c)(1)(ii) of the Rights in Technical Data and Computer Software clause at DFARS 52.227-7013 and/or in similar or successor clauses in the FAR, or the DOD or NASA FAR supplement. Unpublished-rights reserved under the Copyright Laws of the United States. Contractor/manufacturer is Thermo Fisher Scientific Inc., 5225 Verona Road, Madison, WI. 53711-4495.

13. INJUNCTIVE RELIEF. Licensee recognizes and acknowledges that Licensee's breach or violation of any term of this Agreement, including without limitation the use or disclosure of the Confidential Information by the Licensee in a manner inconsistent with the provisions of this Agreement, may cause Thermo Fisher Scientific irreparable damage for which remedies other than injunctive relief may be inadequate, and Licensee agrees that in any request by Thermo Fisher Scientific to a court of competent jurisdiction for injunctive or other equitable relief seeking to restrain such use or disclosure, Licensee will not maintain that such remedy is not appropriate under the circumstances. Licensee further agrees that in the event such equitable relief is granted in the United States, Licensee will not object to courts in other jurisdictions granting provisional remedies enforcing such United States judgments.

14. PARTIAL INVALIDITY; WAIVER. If any provision of this Agreement or the application thereof to any party or circumstances shall be declared void, illegal or unenforceable, the remainder of this Agreement shall be valid and enforceable to the extent permitted by applicable law. In such event the parties shall use their best efforts to replace the invalid or unenforceable provision by a provision that, to the extent permitted by applicable law, achieves the purposes intended under the invalid or unenforceable provision. No failure of either party to exercise any power or right given either party hereunder or to insist upon strict compliance by either party with its obligations hereunder, and no custom or practice of the parties at variance with the terms hereof shall constitute a waiver of either party's right to demand exact compliance with the terms hereof.

15. AUDIT. Thermo Fisher Scientific reserves the right to audit Licensee's use of the Software at the expense of Thermo Fisher Scientific, and Licensee shall allow Thermo Fisher Scientific upon reasonable notice access to Licensee's premises for such purpose. In the event it is determined pursuant to the audit that Licensee's use of the Software is outside the scope of the terms and conditions set forth in this Agreement or Licensee has breached the Agreement, in addition to any other rights or remedies available to Thermo Fisher Scientific hereunder or otherwise, Licensee shall reimburse Thermo Fisher Scientific all the costs of such audit plus any fees owed as a result of the excess use or breach of this Agreement.

16. NOTICES. Any notices to Thermo Fisher Scientific under or in relation to this Agreement shall be sent in writing to: General Manager, Thermo Electron Scientific Instruments LLC, 5225 Verona Road, Madison, Wisconsin 53711-4495.

Detectors for Fourier Transform Spectroscopy

Key Words

- Detector
- FT-IR
- Responsivity
- Sampling Beam Geometry
- Specific Detectivity

Thermo Scientific FT-IR spectrometers can be configured to meet your need to produce the highest performance spectral data for a wide variety of applications. Optimization for each application requires the selection of specific system components – source, optics, electronics, interferometer, beamsplitter, and detector – to produce the best quality spectrum. Of these components, the detector is key, as it can be matched uniquely to specific applications. The choice of the ideal detector for spectral measurement is dependent upon many factors, including:

- Sensitivity
- Spectral range of the measurement
- Sampling beam geometry
- Temporal resolution of the data collection
- Spectral resolution
- Response time

This note will describe the unique features of detectors for vibrational spectroscopy, focusing primarily on Fourier transform infrared (FT-IR). It will help guide the choice for the ideal detector for specific applications.

High throughput (greater than 20% of IR beam reaching the IR detector) and static experiments are generally run utilizing a thermal detector, such as deuterated, L-alanine doped triglycine sulfate (DLaTGS), because it gives full, specific detectivity (D^*) in high-flux environments. From an applications perspective, the DLaTGS detector provides linear response over a wide range of FT-IR throughput, which is beneficial in qualitative and quantitative FT-IR sampling. Thermal detectors are generally less effective for kinetic measurements because their signal is inversely proportional to data collection speed. Low-throughput experiments (less than 20% of IR beam reaching the detector) benefit from the use of a quantum detector, such as the mercuric cadmium telluride (MCT) detector. High MCT sensitivity yields good IR signal in a low-flux measurement. Furthermore, the MCT detector exhibits a relatively constant signal versus data-collection speed and is, therefore, ideal for kinetic measurements. A limitation of MCT detectors is that they lose responsivity and D^* at high throughput. From the applications perspective, this means that in high-throughput applications, we must limit

the beam energy reaching the MCT detector to prevent saturation, where the detector responds with a non-linear signal. Limiting the beam energy can be done by using neutral density filters, or optical screens.

Vibrational spectroscopy detector specifications are listed in Table 1. DLaTGS is a high-performance element which provides higher sensitivity than older technology based on deuterated triglycine sulfate (DTGS) detector elements.

Element size is reported in terms of diameter (mm) for circular detector elements, or area (mm^2), with x, y dimensions for quadrangular shaped elements. Ideally, systems are designed to match the beam diameter and detector element size – not over-filling or under-filling the detector element. A significantly over-filled detector element will not measure the full system flux and, consequentially, the signal-to-noise ratio will be diminished. This must be balanced with the fact that detector noise increases with increased detector area. To prevent saturation in high- D^* detectors, reducing the beam intensity by reducing the beam diameter is not effective. Any area of the detector that is saturated will demonstrate saturation in the resulting spectrum. Figure 1 shows examples of an MCT-A, DLaTGS (with Validation Wheel option) and InGaAs detectors.




Figure 1: Vibrational spectroscopy detectors use a pre-aligned, pinned-in-place design and may be repositioned to other system locations with automatic recognition

Detector Element/Window	Element Size	Spectral Range	Typical D* cm Hz ^{1/2} W ⁻¹	Minimum Responsivity	Operating Temperature	Typical Preamplifier Bandwidth
Mid-infrared						
DLaTGS/KBr	1.3 mm	12500-350 cm ⁻¹	2.7 E8	50 V/W	Room temp.	75 kHz
DLaTGS/KBr	1.3 mm	12500-350 cm ⁻¹	2.7 E8	50 V/W	TE cooled	75 kHz
DLaTGS/CsI	1.3 mm	6400-200 cm ⁻¹	2.4 E8	50 V/W	Room temp.	75 kHz
LiTaO ₃ /KBr	1.5 mm	10000-400 cm ⁻¹	1.0 E8	20 V/W	Room temp.	75 kHz
MCT-A*/CdTe	1.0 x 1.0 mm ²	11700-800 cm ⁻¹	6.4 E10	1,200 V/W	Liquid N ₂	175 kHz
MCT-A/CdTe	1.0 x 1.0 mm ²	11700-600 cm ⁻¹	4.7 E10	750 V/W	Liquid N ₂	175 kHz
MCT-B/KRS-5	1.0 x 1.0 mm ²	11700-400 cm ⁻¹	8.0 E9	50 V/W	Liquid N ₂	175 kHz
Near-infrared						
InGaAs/glass	1.0 mm	12000-3800 cm ⁻¹	4.4 E10	1.3 A/W	Room temp.	170 kHz
InGaAs/glass	1.0 mm	12000-3800 cm ⁻¹	3.0 E11	1.4 A/W	TE cooled	50 kHz
PbSe/sapphire	1.0 x 1.0 mm ²	11000-2000 cm ⁻¹	2.5 E9	6,000 V/W	Room temp.	100 kHz
InSb/CdTe	2.0 x 2.0 mm ²	10000-1850 cm ⁻¹	2.2 E11	2.2 A/W	Liquid N ₂	250 kHz
Visible						
Si/Quartz	2.5 mm	27000-8600 cm ⁻¹	2.8 E12	0.4 A/W	Room temp.	170 kHz
DLaTGS/Quartz	1.3 mm	25000-2000 cm ⁻¹	2.7 E8	50 V/W	Room temp.	170 kHz
Far-infrared						
DLaTGS/Poly	1.5 x 1.5 mm ²	700-50 cm ⁻¹	4.8 E8	300 V/W	Room temp.	15 kHz
Si/Teflon [®]	2.5 x 2.5 mm ²	600-20 cm ⁻¹	1.6 E12	170,000 V/W	Liquid He	1 Hz
IR Microscope						
MCT-A*/CdTe	0.05 x 0.05 mm ²	11700-700 cm ⁻¹	6.4 E10	60,000 V/W	Liquid N ₂	175 kHz
MCT-A*/CdTe	0.25 x 0.25 mm ²	11700-750 cm ⁻¹	6.4 E10	7,000V/W	Liquid N ₂	175 kHz
MCT-A/CdTe	0.25 x 0.25 mm ²	11700-600 cm ⁻¹	4.7 E10	6,500 V/W	Liquid N ₂	175 kHz
MCT-B/KRS-5	0.25 x 0.25 mm ²	11700-450 cm ⁻¹	8.0 E9	400 V/W	Liquid N ₂	175 kHz
InGaAs/glass	0.30 mm	12000-3800 cm ⁻¹	3.0 E12	1.0 A/W	TE cooled	50 kHz
FT-Raman						
InGaAs/glass	1.0 mm	3600-100 R. shift	1.0 E12	.95 A/W	Room temp.	1 kHz
Ge/glass	3.0 mm	3600-100 R. shift	5.0 E13	.70 A/W	Liquid N ₂	2 kHz
PAS						
Photoacoustic	8 mm	4000-100 cm ⁻¹	2.0 E7	5 V/W	Room temp.	30 kHz
TRS						
MCT-A/CdTe	1.0 x 1.0 mm ²	11700-600 cm ⁻¹	4.0 E10	6.0 A/W	Liquid N ₂	20 MHz
MCT-A/CdTe	0.5 x 0.5 mm ²	11700-600 cm ⁻¹	3.4 E10	4.8 A/W	Liquid N ₂	50 MHz

Table 1: Vibrational spectroscopy detector specifications. These specifications are subject to change based upon technological advancement. Please contact your local Thermo Scientific representative for more information.

The spectral ranges for different detectors listed in Table 1 are reported as the endpoints where spectral response begins and ends, and can be diminished based upon specific experimental conditions. Reported values are generated from the combination of detector element responsivity and window material composition. Measured response will be dependent upon the choice of beamsplitter and source and the optics of any sampling accessory used.

Specific D^* is a measure of the detector signal as a function of energy flux and detector noise. At low-light levels, the D^* number may be directly compared from one detector type to another. For example, the MCT-A* detector with a D^* of 6.4×10^{10} is over 200 times more sensitive than the DLaTGS detector with a D^* of 2.7×10^8 . It is important to consider that detectors with high D^* values tend to demonstrate saturation effects in high-throughput experiments. In this case, a neutral-density screen may be used to reduce energy flux at the detector.

D^* is proportional to the signal from the detector relative to the energy flux reaching the detector. This means that a higher D^* detector can provide higher response in a low energy flux environment. D^* accounts for detector noise, so selecting a high D^* detector can provide a higher signal and lower noise.

Figure 2 shows a graphical representation of D^* and the spectral range for the detectors. This figure can be used to graphically compare detector D^* and spectral range for most of the detector products.

Another measure of detector performance is responsivity (R_y), which is a measure of the detector signal per energy flux and detector area. A higher value of R_y indicates a higher detector response for a given energy flux and is corrected for detector area. Note that R_y does not consider detector noise. However, responsivity is a good measure of detector response where signal is sufficient such that noise is not observed.

Many FT-IR detectors operate at room temperature while others require cooling with liquid nitrogen (77 K) for example. Thermo Scientific electric cooling (TE or Peltier cooling) provides enhanced, 100% line stability for the DLaTGS detector, minimizing response changes due to changes in room temperature. For indium gallium arsenide (InGaAs) detectors, cooling can provide increased performance at very low-light levels. TE cooling for DLaTGS improves stability, but not D^* ; cooling quantum detectors improves D^* .

Our detectors cooled by liquid nitrogen utilize a patented, stainless steel dewar, which provides a liquid nitrogen hold time of 18 hours (U.S. Patent 4,740,702). This extended hold time is useful for long experiments such as microscope mapping and extended kinetic measurements. This proprietary detector design also eliminates the formation of ice on the MCT element, thereby preventing spectral artifacts due to water absorption.

Generally, higher bandwidths are required for fast, data-collection applications, such as rapid-scan analysis and time-resolved spectroscopy (TRS).

We offer a wide variety of detectors to optimize each FT-IR experiment and sampling condition. This document provides a reference of general sensitivity, responsivity, and spectral range information for most of the detectors. It clearly demonstrates a wide variety of choices with which you can configure an instrument to precisely meet your needs. It is intended as a general reference for standard versions of the detector products. If you have questions or comments about a detector not shown here, or a special requirement, please contact your local Thermo Scientific representative.




Figure 2: Graphical representation of detector D^* versus spectral range for Nicolet FT-IR detectors. The vertical line at 15798 cm^{-1} represents the helium neon (HeNe) laser emission band which is utilized for internal registration of spectral data.

In addition to these offices, Thermo Fisher Scientific maintains a network of representative organizations throughout the world.

- Africa**
+43 1 333 5034 127
- Australia**
+61 2 8844 9500
- Austria**
+43 1 333 50340
- Belgium**
+32 2 482 30 30
- Canada**
+1 800 530 8447
- China**
+86 10 5850 3588
- Denmark**
+45 70 23 62 60
- Europe-Other**
+43 1 333 5034 127
- France**
+33 1 60 92 48 00
- Germany**
+49 6103 408 1014
- India**
+91 22 6742 9434
- Italy**
+39 02 950 591
- Japan**
+81 45 453 9100
- Latin America**
+1 608 276 5659
- Middle East**
+43 1 333 5034 127
- Netherlands**
+31 76 587 98 88
- South Africa**
+27 11 570 1840
- Spain**
+34 914 845 965
- Sweden/Norway/Finland**
+46 8 556 468 00
- Switzerland**
+41 61 48784 00
- UK**
+44 1442 233555
- USA**
+1 800 532 4752

www.thermo.com



Thermo Electron Scientific Instruments LLC, Madison, WI
USA is ISO Certified.

AN50808_E 11/07M

©2007 Thermo Fisher Scientific Inc. All rights reserved. All trademarks are the property of Thermo Fisher Scientific Inc. and its subsidiaries.

Specifications, terms and pricing are subject to change. Not all products are available in all countries. Please consult your local sales representative for details.

Declaration of Conformity

We,

Thermo Fisher Scientific
Molecular and Microanalysis Products
5225 Verona Rd.
Madison, WI 53711
Ph: +1 608 276 6100

declare under our sole responsibility that the product

Nicolet™ iS™20 FT-IR spectrometer
Nicolet™ iZ10 module

To which this declaration relates is in conformity with the following standard(s) or other normative document(s)

EN 61010-1:2010 (Safety requirements for electrical equipment for measurement, control, and laboratory use)

EN 60825-1:2014 (Safety of laser products)

EN 61326-1:2013 (Electrical equipment for measurement, control, and laboratory use. EMC requirements)

EN 50581:2012 (Assessment of electrical and electronic products with respect to the RoHS)

following the provisions of Directive(s) 2014/30/EU (Electromagnetic Compatibility), 2014/35/EU (Low Voltage), and 2011/65/EU (RoHS).

The Technical Construction file is maintained at:

Thermo Fisher Scientific
Molecular and Microanalysis Products
5225 Verona Rd.
Madison, WI 53711
Ph: +1 608 276 6100

The authorized representative located within the Community is:

Thermo Fisher Scientific
Dieselstr. 4
76227 Karlsruhe, Germany
Ph: +49 (0) 721 4094 389

Date of issue: 04 DEC 2018
 (DD) (MMM) (YYYY)

Place of issue: Madison, Wisconsin, USA

Signed: _____


Carlo Nardini
Director of Operations



ThermoFisher
SCIENTIFIC